

# T4 DNA Polymerase

## (5,000 U/ml)

Error Rate:  $\sim 1 \times 10^{-6}$  bases



**Catalog:** RK20539

**Size:** 200 U / 1000 U

**Concentration:** 5,000 U/ml

**Components:**

T4 DNA Polymerase (5,000 U/ml)	RM21302
10X ABuffer B	RM20126

## Product Description

T4 DNA Polymerase catalyzes the synthesis of DNA in the 5' → 3' direction and requires the presence of template and primer. This enzyme has a 3' → 5' exonuclease activity, which is much more active than that found in DNA Polymerase I (*E. coli*). Unlike *E. coli* DNA Polymerase I, T4 DNA Polymerase does not have a 5' → 3' exonuclease function.

It is applicable to 3' overhang removal to form blunt ends, 5' overhang fill-in to form blunt ends, single-strand deletion subcloning, second-strand synthesis in site-directed mutagenesis, and probe labeling using replacement synthesis.

**Product Source:** Purified from a strain of *E. coli* that carries the T4 DNA Polymerase gene.

**Unit Definition:** One unit is defined as the amount of enzyme that incorporates 10 nmol of dNTP into acid insoluble material in 30 minutes at 37°C.

**Storage Conditions:** 100 mM KPO<sub>4</sub>, 1 mM DTT, 50% Glycerol, pH 6.5 @ 25°C

**Storage Temperature:** -20°C

**Reaction Conditions:** 1X ABuffer B, Incubate at 12°C.

**1X ABuffer B:** 10 mM Tris-HCl, 50 mM NaCl, 10 mM MgCl<sub>2</sub>, 1 mM DTT, pH7.9 @ 25°C

**Activity in ABclonal Buffer**

ABufferA	ABufferB	ABufferC	ABufferS
60%	100%	100%	100%

**Heat Inactivation:** 75°C for 20 min

**Molecular Weight:** Theoretical 104000 daltons

**5' - 3' Exonuclease:** No

**3' - 5' Exonuclease:** Yes

**Strand Displacement:** No

**Order:** [service@abclonal.com](mailto:service@abclonal.com)  
**Tech:** [support@abclonal.com](mailto:support@abclonal.com)

For research purposes only. Not for therapeutic or diagnostic purposes in humans or animals.  
Please visit <http://abclonal.com> for a complete listing of recommended products.

## Instructions

**Protocol for blunting ends by 3' overhang removal and 3' recessed (5' overhang) end fill-in using T4 DNA Polymerase:**

1. DNA should be dissolved in 1X reaction buffer\* supplemented with 100 μM dNTPs.
2. Add 1 unit T4 DNA Polymerase per microgram DNA.
3. Incubate 15 minutes at 12°C.
4. Stop reaction by adding EDTA to a final concentration of 10 mM and heating to 75°C for 20 minutes.

**CAUTION:** Elevated temperatures, excessive amounts of enzyme, failure to supplement with dNTPs, or long reaction times will result in recessed ends due to the 3' → 5' exonuclease activity of the enzyme.

\*T4 DNA Polymerase can be used in *CutA*, *CutB*, and *CutS* as well as ABuffer A/B/S and T4 DNA Ligase Reaction Buffer. Optimal activity is observed in ABufferB. BSA supplementation is recommended when using a buffer that does not already contain BSA.

### Notes:

1. For fill-in reactions only: T4 DNA Polymerase can be used in *CutA*, *CutB*, and *CutS* as well as ABuffer A/B/S and T4 DNA Ligase Reaction Buffer.
2. For blunting reactions requiring removal of overhangs: T4 DNA Polymerase can be used in *CutA*, *CutB*, and *CutS* as well as ABuffer A/B/S and T4 DNA Ligase Reaction Buffer. *CutC* and ABuffer C are not recommended when overhang removal is required.
3. Optimal activity is observed in ABuffer B.
4. Supplement with dNTPs\*.
5. BSA supplementation is recommended when using a buffer that does not already contain BSA.
6. Incubate at temperature suggested for specific protocol.

\* Refer to specific protocol to determine recommended dNTP concentrations.

### QC Process:

- Purity is above 95% detected by SDS-PAGE.
- No nuclease, RNase contamination.
- No residual host genomic DNA is detected by PCR.