

LongerAmp Taq DNA Polymerase



Catalog: RK20680

Size: 1,000 U / 2,500U

Concentration: 2,500 U/ml

Components:

LongerAmp Taq DNA Polymerase (2,500 U/ml)	RM20312
5X LongerAmp Taq Reaction Buffer	RM20140

Product Description

LongerAmp Taq DNA Polymerase is a unique blend of Taq klenow and Deep Ocean DNA Polymerases. The 3'→5' exonuclease activity of Deep Ocean DNA Polymerase increases the fidelity and robust amplification of Taq klenow DNA Polymerase. LongerAmp Taq DNA Polymerase offers two fold higher fidelity than Taq klenow DNA Polymerase alone. A wide range of PCR products can be generated, up to 30 kb from lambda DNA or from human genomic DNA.

Unit Definition:

One unit is defined as the amount of enzyme that will incorporate 10 nmol of dNTP into acid insoluble material in 30 minutes at 75°C.

Storage Conditions:

10 mM Tris-HCl, 100 mM KCl, 1 mM DTT, 0.1 mM EDTA, 0.5% Tween 20, 0.5% IGEPAL CA-630, 50% Glycerol, pH7.4 @ 25°C

Storage Temperature: -20°C

Reaction Conditions: 1X LongerAmp Taq Reaction Buffer

1X LongerAmp Taq Reaction Buffer:

60 mM Tris-SO₄, 20 mM (NH₄)₂SO₄, 2 mM MgSO₄, 3% Glycerol, 0.06% IGEPAL CA-630, 0.05% Tween 20, pH9.1 @ 25°C

Application Features:

Long Range PCR; Colony PCR

Heat Inactivation: No

5' - 3' Exonuclease: No

3' - 5' Exonuclease: Yes

Strand Displacement: +

Resulting Ends: Mix of Blunt and Single-base 3' Overhangs

Notes: 5'→3' flap endonuclease destroys displaced strand.

Instructions

Reaction setup: We recommend assembling all reaction components on ice and quickly transferring the reactions to a thermocycler preheated to the denaturation temperature (94°C).

Component	25 µl	50 µl	Final Conc.
Nuclease-free water	to 25 µl	to 50 µl	
5X LongerAmp Taq Reaction Buffer	5 µl	10 µl	1X
10 mM dNTPs	0.75 µl	1.5 µl	300 µM
10 µM Forward Primer	1 µl	2 µl	0.4 µM (0.05–1 µM)
10 µM Reverse Primer	1 µl	2 µl	0.4 µM (0.05–1 µM)
Template DNA	variable	variable	<1,000 ng
LongerAmp Taq DNA Polymerase	1 µl	2 µl	5 units/50 µl PCR

Notes: Gently mix the reaction. Avoid pipetting samples containing target DNA when amplicons above 20 kb are desired. Collect all liquid to the bottom of the tube by a quick spin if necessary. Overlay the sample with mineral oil if using a PCR machine without a heated lid.

Transfer PCR tubes from ice to a PCR machine with the block preheated to 94°C and begin thermocycling. A PCR program was recommended.

Thermocycling conditions for a routine PCR:

STEP	STEP	TIME
Initial Denaturation	94°C	3 minutes
30 Cycles	98°C	5-10seconds
	55-60°C(optional)	20-30seconds
	65-68°C	50 seconds/kb
Final Extension	65-68°C	10 minutes
Hold	4-10°C	

Order: order@abclonal.com
Tech: support@abclonal.com

For research purposes only. Not for therapeutic or diagnostic purposes in humans or animals.
Please visit <http://abclonal.com.cn> for a complete listing of recommended products.

General Guidelines:

1. Template:

The quality of the DNA template is essential for long-range PCR amplification. Recommended amounts of DNA template for a 50 µl reaction are as follows:

DNA	up to 15 kb	above 15 kb
genomic	1 ng–500 ng	10 ng–1 µg
plasmid or viral	1 pg–1 ng	10 pg–10 ng

Successful amplification above 20 kb largely depends on the quality of DNA templates and the primer sequences.

2. Primers:

Oligonucleotide primers are generally 20–40 nucleotides in length and ideally have a GC content of 40–60%. For amplicons larger than 20 kb, it is desirable to have primers with GC content above 50%, matched T_m above 60°C and primers at least 24 nucleotides in length. The final concentration of each primer in a PCR reaction may be 0.05–1 µM, typically 0.1–0.5 µM.

3. Mg⁺⁺ and additives:

Mg⁺⁺ concentration of 1.5–2.0 mM is optimal for most PCR products generated with LongerAmp Taq DNA Polymerase. The final Mg⁺⁺ concentration in 1X LongerAmp Taq Reaction Buffer is 2 mM. This supports satisfactory amplification of most amplicons. However, Mg⁺⁺ can be further optimized in 0.5 or 1.0 mM increments using MgSO₄.

Amplification of some difficult targets, like GC-rich sequences, may be improved with additives, such as DMSO or formamide.

4. Deoxynucleotides:

The recommended final concentration of dNTPs for long-range PCR is 300 µM of each deoxynucleotide.

5. LongerAmp Taq DNA Polymerase concentration:

We generally recommend using LongerAmp Taq DNA Polymerase at a concentration of 100 units/ml (5 units/50 µl reaction). However, the optimal concentration of LongerAmp Taq DNA Polymerase may vary in specialized applications.

6. Denaturation:

An initial denaturation of 3 minutes at 94°C is sufficient for most amplicons from pure DNA templates. For difficult templates or longer amplicon, a higher initial denaturation at 98°C is recommended prior to PCR cycling to fully denature

the template.

During thermocycling a 5–10 second denaturation at 98°C is recommended.

7. Annealing:

The annealing step is typically 15–60 seconds. Annealing temperature is based on the T_m of the primer pair and is typically 45–65°C. Annealing temperatures can be optimized by doing a temperature gradient PCR starting 5°C below the calculated T_m. When primers with annealing temperatures above 60°C are used, a 2-step PCR protocol is possible (see #10).

8. Extension:

The recommended extension temperature is 65–68°C. Extension times are generally 50 seconds per kb. A final extension of 10 minutes at 65–68°C is recommended. Notice: 72°C was not suitable for this product.

9. Cycle number:

Generally, 25–35 cycles yields sufficient product. Up to 45 cycles may be required to detect low-copy-number targets.

10. 2-step PCR:

When primers with annealing temperatures above 60°C are used, a 2-step thermocycling protocol is possible. Moreover, When PCR product is longer than 5kb, a 2-step PCR protocol is recommended especially for human or *E.coli* genomic.

Thermocycling conditions for a routine 2-step PCR:

Temperature	Time	Cycles
94°C	3min	1
98°C	5-10s	30 cycles
65-68°C	50s/kb	
65-68°C	10 min	1
4-10°C	∞	

11. PCR product:

The majority of the PCR products generated using LongerAmp Taq DNA Polymerase contain dA overhangs at the 3'-end; therefore the PCR products can be ligated to dT/dU-overhang vectors.

12. Optimization:

Denature temperature, extension temperature, and different time for every step were optimized for special and difficult amplicon. This optimization was very important for successful results.

Order: order@abclonal.com *For research purposes only. Not for therapeutic or diagnostic purposes in humans or animals.*
Tech: support@abclonal.com Please visit <http://abclonal.com.cn> for a complete listing of recommended products.