




# Rapid Plus DNA Lib Prep kit

## RK20208

---

(1ng-1 $\mu$ g Input DNA)

(Illumina<sup>®</sup> Compatible)



[www.abclonal.com](http://www.abclonal.com)  
version: N12G13v1.0

# Contents

1.Introduction .....	01
2.List of Components .....	02
3.Additional Materials Required .....	02
4.Precautions .....	03
5.Protocol .....	03
Step 1. End Preparation .....	03
Step 2. Adaptor Ligation .....	04
Step 3. Size Selection (optional) .....	05
Step 4. Amplification (optional) .....	06
6.Primer Sequences .....	08


---

# 1.Introduction

The Rapid Plus DNA Lib Prep kit provides a streamlined and efficient method to construction DNA library for next generation sequencing (NGS) at Illumina® platforms. The entire protocol can be completed in 2 hours. Starting from fragmented input-DNA (1 - 1000 ng), the kit contains all the enzymes and buffers for end-preparation, adaptor ligation and library amplification. To simplify the workflow, the fragmented DNA is end-repaired and adenylated during the initial step. The kit can be conveniently used for construction of PCR-free library from 100 ng of high-quality fragmented genomic DNA, or 10 ng of circulating cell-free DNA (cfDNA). PCR primers provided in the kit are optional for users who prefer self-selected high-fidelity DNA polymerases in the library amplification.

The kit also provides DNA repair enzymes, which could repair DNA damages such as nicks, gaps, oxidized bases, damaged/blocked 3' ends, AP sites (apurinic/apyrimidinic sites), and uracil bases. In addition, the DNA repair enzymes significantly increase the NGS library yields from FFPE-derived DNA.

The following workflow illustrates the processes of Rapid DNA Lib Prep Kit.



## 2. List of Components

All components should be stored at -20°C. The shelf lives of all reagents are one year when stored properly.

**Table 1. Kit Contents**

	Tube name	8 RXN	24 RXN	96 RXN
End Preparation	End Prep Buffer	70 µl	168 µl	700 µl
	End Prep Enzymes	30 µl	78 µl	300 µl
	Damaged DNA Repair Enzymes	30 µl	78 µl	300 µl
Adaptor Ligation	Ligation Buffer	300 µl	780 µl	3 ml
	Ligation Enzymes	100 µl	240 µl	1000 µl
Amplification	2X PCR Master Mix	250 µl	650 µl	2.5 mL
	10X PCR Primers	50 µl	130 µl	500 µl

## 3. Additional Materials Required

Mechanical or enzymatic methods for DNA fragmentation

Multiplex adapters (cat. no. RK20282, RK20283, RK20284, RK20285, RK20292, RK20293, RK20294, RK20295) compatible with Illumina® platforms

100% ethanol (ACS grade)

Nuclease-free water

PCR strip tubes or plates

Magnetic stand

Thermocycler

Agencourt™ AMPure XP bead (Beckman Coulter Inc., cat. no. A63880)

Pipettes and multichannel pipettes

Aerosol resistant pipette tips

Microcentrifuge

Vortex mixer

Agilent Bioanalyzer or comparable method to assess the quality of DNA library

## 4. Precautions

- Make sure to use high-quality DNA samples. Heavily nicked or damaged DNA will significantly lower library yield. We recommend Qubit® or other fluorometric methods to quantify DNA input. Residual trace RNA, nucleotides, ss DNA or other contaminants will decrease the efficiency of enzymes for optimized library preparation.
- Do not vortex enzymes, mix components by gently pipetting up and down several times.

## 5. Protocol

- Size selection step is optional for DNA libraries obtained after adapter ligation. For DNA input less than 50 ng, the size selection is not recommended.
- PCR-based amplification step is normally required if the DNA input amount is less than 50 ng for genomic DNA and 10 ng for cfDNA.
- For DNA input less than 10 ng, an additional cleanup step is highly recommended after ligation reaction.
- The full length Y-shape adapters (cat. no. RK20282, RK20283, RK20292, RK20293) suitable for Illumina® sequencing platforms are sold separately.
- If using high-fidelity DNA polymerase from other suppliers, PCR primers are included in the kit.

### Step 1. End Preparation

- **Prepare fresh 80% ethanol**

1.1 Prepare end-preparation reaction mix in PCR tubes on ice according to the Table 2. below.

**Table 2. End-preparation Reaction Setup (per sample)**

Fragmented DNA	X µl
End Prep Buffer*	7 µl
End Prep Enzymes*	3 µl
Damaged DNA Repair Enzymes*	3 µl
Nuclease-free Water	Up to 60 µl
Total Volume	60 µl

\* The buffer and enzyme mix should be pre-mixed and then added in a single pipetting step.

1.2 Mix thoroughly by pipetting.

1.3 Incubate reaction tubes in a thermocycler according to the program listed in Table 3. (**WITH** heated lid set at 75°C).

**Table 3. End-preparation Reaction Program**

30°C	30 min
65°C	30 min
4°C	∞

## Step 2. Adaptor Ligation

2.1 Prepare and dilute adapters in low-EDTA TE buffer or nuclease-free water according to the instructions in the Table 4. Below.

**Table 4. Adaptor Dilution**

Input DNA	Adaptor Dilution	Adapter Concentration
1 µg-50 ng	No dilution	15 µM
49 ng-25 ng	2-fold	7.5 µM
24 ng-10 ng	5-fold	3 µM
9 ng-5 ng	10-fold	1.5 µM
<5 ng	20-fold	0.75 µM

2.2 Prepare the ligation reaction mix in PCR tubes on ice according to the Table 5. Below.

**Table 5. Ligation Reaction Setup**

End Prep Reaction Mix (Step1.3)	60 µl
H <sub>2</sub> O*	5 µl
Ligation buffer*	30 µl
Working Adaptor	5 µl
Ligation Enzymes*	10 µl
Total Volume	110 µl

\* The water, buffer and ligation enzymes should be premixed and then added in a single pipetting step.

- 2.3 Incubate at 20°C for 15 minutes in a thermocycler **WITHOUT** a heated lid, and then hold at 4°C.
- 2.4 Clean up ligated DNA.
  - 2.4.1 Add 88 µl (ratio 0.8X) of Agencourt™ Ampure XP beads to each samples, and mix well by pipetting.
  - 2.4.2 Incubate the mixture at room temperature (RT) for 5 minutes.
  - 2.4.3 Pellet the beads on a magnetic stand at RT for 2 minutes
  - 2.4.4 Carefully remove and discard the supernatant.
  - 2.4.5 Wash the beads with 200 µl fresh 80% ethanol. Pellet the beads on a magnetic stand and carefully remove the ethanol.
  - 2.4.6 Repeat step 2.4.5 for a total of two washes
  - 2.4.7 Resuspend the magnetic beads in 21 µl nuclease-free water for **amplification**. Or resuspend the magnetic beads in 51 µl nuclease-free water for **double-size selection**.
  - 2.4.8 Mix thoroughly by pipetting, and then incubate at RT for 1 minute to release the DNA from the beads.
  - 2.4.9 Pellet the beads on a magnetic stand at RT for 2 minutes.
  - 2.4.10 Transfer 20 µl or 50 µl of the supernatant to a new PCR tube.
  - 2.4.11 Store the library at -20°C until ready for QC, library quantification or sequencing.

### Step 3. Size Selection (optional)

**Table 6. Ratios of Agencourt™ AMPure XP Beads for DNA Size Selection**

Median Insert Size(bp)	150~350
Library Size (bp)	250~500
1 <sup>st</sup> Binding Beads	0.7x (35 µl)
2 <sup>nd</sup> Binding Beads	0.2x (10 µl)

- 3.1 Add 1<sup>st</sup> binding beads into 50 µl of adapter ligated DNA (from step 2.4.10), according to the volume ratio described in the Table 6., and mix thoroughly by pipetting.
- 3.2 Incubate at RT for 5 minutes.
- 3.3 Pellet the beads on a magnetic stand at RT for 2 minutes.
- 3.4 Carefully transfer the supernatant to a new PCR tube.
- 3.5 Add 2<sup>nd</sup> binding beads to the supernatant, according to the volume ratio described in the Table 6., and mix thoroughly by pipetting.
- 3.6 Incubate at RT for 5 minutes
- 3.7 Pellet the beads on a magnetic stand at RT for 2 minutes.
- 3.8 Carefully remove and discard the supernatant.
- 3.9 Wash the beads with 200 µl fresh 80% ethanol. Pellet the beads on a magnetic stand and carefully remove the ethanol.
- 3.10 Repeat Step 3.9 for a total of two washes.
- 3.11 Air dry the beads for 5 minutes on a magnetic stand.
- 3.12 Resuspend magnetic beads in 21 µl nuclease-free water.
- 3.13 Mix thoroughly by pipetting, and then incubate at RT for 1 minute to release the DNA from the beads.
- 3.14 Pellet the beads on a magnetic stand at RT for 2 minutes.
- 3.15 Transfer 20 µl of the supernatant to a new PCR tube for amplification.

## **Step 4. Amplification (optional)**

- 4.1 Prepare the PCR reaction according to the Table 7. below.

**Table 7. PCR Amplification Reaction Setup**



Adapter-Ligated DNA	20 $\mu$ l
2x PCR Master Mix	25 $\mu$ l
10x PCR Primers	5 $\mu$ l
<b>Total Volume</b>	<b>50 <math>\mu</math>l</b>

4.2 Mix thoroughly by pipetting.

4.3 Program a thermocycler according to the Table 8. below.

**Table 8. PCR Cycles for Library Amplification**

98°C	1 min	
98°C	10 s	
60°C	30 s	2-17 PCR cycles
72°C	30 s	
72°C	1 min	
12°C	$\infty$	

4.4 Add 50  $\mu$ l (ratio 1.0X) of Agencourt™ AMPure XP beads to each reaction tubes, mix thoroughly by pipetting.

4.5 Incubate at RT for 5 minutes.

4.6 Pellet the beads on a magnetic stand at RT for 2 minutes.

4.7 Carefully remove and discard the supernatant.

4.8 Wash the beads with 200  $\mu$ l fresh 80% ethanol. Pellet the beads on a magnetic stand and carefully remove the ethanol.

4.9 Repeat Step 4.8 for a total of two washes.

4.10 Resuspend the magnetic beads in 21  $\mu$ l of low-EDTA TE buffer.

4.11 Mix thoroughly by pipetting, and then incubate at RT for 1 minute to release the DNA from the beads.

4.12 Pellet the beads on a magnetic stand at RT for 2 minutes.

4.13 Transfer 20  $\mu$ l of clear supernatant to a new PCR tube.

4.14 Store the library at -20°C until ready for QC, library quantification or sequencing.

## 6. Primer Sequences

### Primer 1:

5'-AATGATACGGCGACCACCGAG

### Primer 2:

5'-CAAGCAGAAGACGGCATACGAG



**United States**

[www.abclonal.com](http://www.abclonal.com)

Address: 86 Cummings Park Dr ,Woburn,MA 01801, United States

Phone: 888.754.5670, +1 857.259.4898 (Int'l)

Email: [service@abclonal.com](mailto:service@abclonal.com)