

Immunohistochemistry on paraffin-embedded Protocol

- 1.** Deparaffinizing twice, 10 min each time.
- 2.** Rehydration
- 3.** Antigen retrieval and then natural cooling
- 4.** Primary antibody incubation. Incubate slides with 3 % H₂O₂ solution (100μL) for 15 min to quench endogenous peroxidase activity. Rinse slides with PBS, 3 times and 5 min each time.
- 5.** Draw 2 lines beside tissue section to keep the liquid from flowing away.
- 6.** Block the sections with corresponding serum or BSA.
- 7.** Incubate sections with primary antibody overnight at 4 °C. Then rinse slides with PBS, 3 times and 5 min each time.
- 8.** Apply peroxidase labeled secondary antibody and incubate for 30 min. Rinse slides with PBS.
- 9.** DAB plus chromogen
- 10.** Hematoxylin or DAPI for nuclei staining.
- 11.** Dehydration and mounting. Dehydrate sections with graded ethanol, 75 % (1 min), 85 % (1 min), 90 % (1 min), 100 % (1 min), 100 % (1 min). Immerse slides in xylene for 5 min. Repeat once.
- 12.** Signal Detection. Mount and cover the section for detection.

Immunohistochemistry - Frozen Protocol

1. Wash sections in 1X TBS two times for 5 min.
2. Incubate for 10 min at room temperature in methanol/peroxidase.
3. Wash sections in 1X TBS two times for 5 min.
4. Block each section with 100–400 μ l 1X TBS/0.3% Triton X-100/5% Normal Goat Serum for 1 hour at room temperature.
5. Remove blocking solution and add 100–400 μ l primary antibody diluted in 1X TBS/0.3% Triton X-100/5% Normal Goat Serum to each section.
6. Incubate overnight at 4°C.
7. Return to room temperature.
8. Remove antibody solution and wash sections in 1X TBS three times for 5 min each.
9. Incubate in a humidified chamber for 30 min at room temperature.
10. Wash sections three times with 1X TBS for 5 min each.
11. Prepare an appropriate volume of substrate solution prior to use by mixing one drop of Liquid DAB plus chromogen immediately with 1 ml of substrate buffer. Apply the substrate carefully and incubate for 5-10 minutes till a brown color develops.
12. Immerse slides in dH₂O.
13. To stain nuclei, immerse slides in a bath of hematoxylin for 3 minutes.
14. Rinse slides gently with distilled water.
15. Transfer slides into a 1% HCl, 99% ethanol solution for 10 seconds; transfer to distilled water immediately.
16. Wash sections in dH₂O two times for 5 min each.
17. Dehydrate sections: Incubate sections in 95% ethanol two times for 10 sec each. Repeat in 100% ethanol, incubating sections two times for 10 sec each. Repeat in xylene, incubating sections two times for 10 sec each.
18. Mount sections with coverslips.