




2X Frag/Elute Buffer[®]

RM20225



www.abclonal.com

Version: N17 H14v1.1

Table of Contents

| | |
|---------------------|---|
| 1. Product Overview | 1 |
| 2. Kit Components | 2 |
| 3. Storage | 2 |
| 4. Protocol | 4 |

1. Product Overview

- ◆ 2X Frag/Elute Buffer (ABclonal, Cat. RM20225) was used for RNA Fragmentation and First-strand synthesis steps in RNA library construction.
- ◆ The kit is applicable to total RNA samples or purified mRNA samples from eukaryotes, including animals, plants, and fungi.
- ◆ 2X Frag/Elute Buffer contains Mg^{2+} and interrupts the RNA by high temperature, so the RNA samples used for fragmentation cannot contain metal ion components such as Mg^{2+} and EDTA.
- ◆ 2X Frag/Elute Buffer contains Random Primer as reverse transcription primer in the First strand cDNA synthesis step.

2. Kit Components

| Tube Name | 96 RXN |
|----------------------|--------|
| 2X Frag/Elute Buffer | 528 µL |

3. Storage

- ◆ Storage: -20°C to -10°C.
- ◆ Long-distance transportation: The kit should be transported with dry ice or with both dry ice and ice bags at -40°C to -20°C

4. Protocol

1. RNA Fragmentation

1.1 Prepare the Prepare the RNA fragmentation system according to the following table:

| Reagent | Volume |
|----------------------|------------------|
| Input RNA | 10-1000 ng |
| 2X Frag/Elute Buffer | 6 μ L |
| Nuclease-free Water | Up to 12 μ L |

1.2 If the Input RNA does not need to be fragmented, 2. First strand cDNA synthesis was performed directly; If the Input RNA needs to be fragmented, please perform the RNA fragmentation procedure according to the recommended conditions in the table below (heating lid temperature 105°C):

| Target Fragment Size | Fragmentation Condition |
|----------------------|-------------------------|
| 200-300 nt | 94°C 15 min, 4°C hold |
| 300-450 nt | 94°C 10 min, 4°C hold |
| 400-700 nt | 94°C 5 min, 4°C hold |

1.3 After cooling to 4°C, take out the tube, centrifuge it instantaneously, and place the tube on the magnetic rack until the solution becomes clear. Transfer 10 μ L of the supernatant into another PCR tube, and immediately use it for the first strand cDNA synthesis.

2. First Strand cDNA Synthesis

2.1 Thaw RT Reagent at room temperature, and prepare the following system on ice:

| Reagent | Volume |
|------------------------------------|------------|
| Fragmented RNAs | 10 μ L |
| RT Reagent* | 8 μ L |
| First Strand Synthesis Enzyme Mix* | 2 μ L |
| Total volume | 20 μ L |

* : The pre-mix can be prepared in advance, and its volume should be 1.1 times the sample volume to make up for the natural loss.

** : The 2X Frag/Elute Buffer contains the Random Primer required for the First strand cDNA synthesis. Ensure that the 2X Frag/Elute Buffer is added when performing this step.

2.2 Mix the prepared system well by pipetting, centrifuge it instantaneously, and incubate it in the PCR system (heating lid temperature 105°C).

| Temperature | Time |
|-------------|--------|
| 25°C | 10 min |
| 42°C | 15 min |
| 70°C | 15 min |
| 4°C | Hold |

China

www.abclonal.com.cn

Headquarters: Building 5, Precision Medicine Industry Base Project I, Gaokeyuan 3rd Road, Donghu New Technology Development Zone, Jiangxia District, Wuhan, Hubei, China

Shanghai R&D Center: F4, Building 2, Zizhu High-Tech Industrial Development Park, No. 58 Yuanmei Road, Minhang District, Shanghai, China

US R&D Center: 86 Cummings Park Dr, Woburn, MA 01801, United States

Tel: 400-999-6126

Email: cn.market@abclonal.com

