

Klenow Fragment (3'→5' exo⁻)



Catalog: RK20526

Size: 200 U/1,000 U/5,000 U

Concentration: 5,000 U/ml

Components:

Klenow Fragment 3'→5' exo ⁻	RM20516
10X ABuffer B	RM20126

Product Description

Klenow Fragment (3' → 5' exo⁻) is an N-terminal truncation of DNA Polymerase I which retains polymerase activity, but has lost the 5' → 3' exonuclease activity and has mutations (D355A, E357A) which abolish the 3' → 5' exonuclease activity.

Klenow Fragment (3' → 5' exo⁻) is isolated from a recombinant source. It generates probes using random primers and shows moderate strand displacement activity.

It is applicable to random priming labeling, DNA sequencing by the Sanger dideoxy method, second strand cDNA synthesis, and second strand synthesis in mutagenesis protocols.

Product Source: An *E.coli* strain containing a plasmid with a fragment of the *E.coli polA* (D355A, E357A) gene starting at codon 324.

Unit Definition: One unit is defined as the amount of enzyme that will incorporate 10 nmol of dNTPs into acid insoluble material in 30 minutes at 37°C.

Storage Conditions: 25 mM Tris-HCl, 1 mM DTT, 0.1 mM EDTA, 50% Glycerol, pH 7.4 @25°C

Storage Temperature: -20°C

Reaction Conditions: 1X ABuffer B, Incubate at 37°C

1X ABuffer B: 10 mM Tris-HCl, 50 mM NaCl, 10 mM MgCl₂, 1 mM DTT, pH7.9 @ 25°C

Heat Inactivation: 75°C for 20 min

Molecular Weight: Theoretical 68000 daltons

5' - 3' Exonuclease: No

3' - 5' Exonuclease: No

Strand Displacement: +++

Error Rate: ~ 100 x 10⁻⁶ bases

Phone: 888.754.5670

Email: service@abclonal.com

*For research purposes only. Not for therapeutic or diagnostic purposes in humans or animals.
For more information, visit www.abclonal.com*

Instructions

A-Tailing with Klenow Fragment (3'→5' exo⁻)

Starting Material:

1-5 µg of blunt-ended DNA* (100-1000 bp).

**If starting with blunt-ended DNA that has been prepared by PCR or by end polishing, DNA must be purified to remove the blunting enzymes.*

1. Mix the following components in a sterile microfuge tube:

Composition	Amount
ddH ₂ O	to 50 µl
10X ABuffer B	5 µl
Purified Blunt DNA	1-5 µg
dATP (10 mM)	0.5 µl (0.1 mM final)
Klenow Fragment (3'→5' exo ⁻)	3 µl (15 U)

2. Incubate in a thermal cycler for 30 minutes at 37°C.
3. Purify DNA sample on one spin column.

Notes:

- Klenow Fragment (3'→5' exo⁻) is not suitable for generating blunt ends because it lacks the 3'→5' exonuclease necessary to remove non-templated 3' additions.
- Klenow Fragment (3'→5' exo⁻) is also active in all ABuffer A/B/C/S when supplemented with dNTPs.
- When Klenow Fragment (3'→5' exo⁻) is used to sequence DNA using the dideoxy method of Sanger *et al.*, 1 unit/5 µl reaction volume is recommended.

QC Process:

- Purity (SDS-PAGE) is above 95%.
- No exonuclease, nuclease, or RNase contamination.
- No residual host genomic DNA detected by PCR.