

A0307

Leader in Biomolecular Solutions for Life Science



# ICAM1 Mouse mAb

Catalog No.: A0307

## Basic Information

### Observed MW

60kDa

### Calculated MW

58kDa

### Category

SMab Recombinant Monoclonal Antibody

### Applications

WB,IF/ICC,ELISA

### Cross-Reactivity

Human

## Background

This gene encodes an integral membrane protein that binds leukocyte adhesion protein LFA-1. It participates in the innate immune response.

## Recommended Dilutions

WB 1:500 - 1:1000

IF/ICC 1:20 - 1:50

## Immunogen Information

### Gene ID

15894

### Swiss Prot

### Immunogen

Recombinant protein of mouse ICAM1

### Synonyms

CD54; Ly-47; Icam-1; MALA-2; ICAM1

## Contact



[www.abclonal.com](http://www.abclonal.com)

## Product Information

### Source

Mouse

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.