

A0597

Leader in Biomolecular Solutions for Life Science



RAB11A Mouse mAb

Catalog No.: A0597

Basic Information

Observed MW

26kDa

Calculated MW

24kDa

Category

SMab Recombinant Monoclonal Antibody

Applications

WB, IHC-P, ELISA

Cross-Reactivity

Human, Mouse, Rat

Background

The protein encoded by this gene belongs to the Rab family of the small GTPase superfamily. It is associated with both constitutive and regulated secretory pathways, and may be involved in protein transport. Two transcript variants encoding different isoforms have been found for this gene.

Recommended Dilutions

WB	1:500 - 1:1000
IHC-P	1:50 - 1:100

Immunogen Information

Gene ID

8766

Swiss Prot

P62491

Immunogen

A synthetic peptide of human RAB11A

Synonyms

YL8; RAB11A

Contact



www.abclonal.com

Product Information

Source

Mouse

Isotype

IgG

Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, pH 7.3.