

A10044

Leader in Biomolecular Solutions for Life Science



ZC3H11A Rabbit pAb

Catalog No.: A10044

Basic Information

Observed MW

120kDa

Calculated MW

89kDa

Category

Mouse Monoclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Enables RNA binding activity. Involved in poly(A)+ mRNA export from nucleus. Colocalizes with transcription export complex.

Recommended Dilutions

WB 1:1000 - 1:4000

Immunogen Information

Gene ID

9877

Swiss Prot

O75152

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human ZC3H11A (NP_055642.3).

Synonyms

ZC3HDC11A; ZC3H11A

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of various lysates using ZC3H11A Rabbit pAb (A10044) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 10s.