

A10362

Leader in Biomolecular Solutions for Life Science



RBM4 Rabbit pAb

Catalog No.: A10362

Basic Information

Observed MW

Refer to figures

Calculated MW

40kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse, Rat

Background

Enables RNA binding activity and cyclin binding activity. Involved in several processes, including IRES-dependent translational initiation of linear mRNA; negative regulation of translation; and regulation of alternative mRNA splicing, via spliceosome. Located in cytoplasmic stress granule; cytosol; and nuclear lumen.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

5936

Swiss Prot

Q9BWF3

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human RBM4 (NP_002887.2).

Synonyms

LARK; RBM4A; ZCRB3A; ZCCHC21; RBM4

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.