

A12978

Leader in Biomolecular Solutions for Life Science



IRF2 Rabbit mAb

Catalog No.: A12978

Recombinant

Basic Information

Observed MW

Refer to Figures

Calculated MW

39kDa

Category

SMab Recombinant Monoclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

IRF2 encodes interferon regulatory factor 2, a member of the interferon regulatory transcription factor (IRF) family. IRF2 competitively inhibits the IRF1-mediated transcriptional activation of interferons alpha and beta, and presumably other genes that employ IRF1 for transcription activation. However, IRF2 also functions as a transcriptional activator of histone H4.

Recommended Dilutions

WB 1:1000 - 1:2000

Immunogen Information

Gene ID

3660

Swiss Prot

P14316

Immunogen

A synthetic peptide of human IRF2

Synonyms

IRF-2; IRF2

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.