

A14867

Leader in Biomolecular Solutions for Life Science



LRP10 Rabbit pAb

Catalog No.: A14867

Basic Information

Observed MW

Refer to figures

Calculated MW

76kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse

Background

This gene encodes a low density lipoprotein receptor family protein. A similar protein in mouse is thought to play a role in the uptake of apolipoprotein E-containing lipoproteins.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

26020

Swiss Prot

Q7Z4F1

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 17-120 of human LRP10 (NP_054764.2).

Synonyms

LRP9; LRP-10; MST087; MSTP087; LRP10

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.