

A16589

Leader in Biomolecular Solutions for Life Science



## DHX57 Rabbit pAb

Catalog No.: A16589

### Basic Information

**Observed MW**

150-155kDa

**Calculated MW**

156kDa

**Category**

Polyclonal Antibody

**Applications**

WB, ELISA

**Cross-Reactivity**

Human, Rat

### Background

Enables RNA binding activity.

### Recommended Dilutions

WB 1:500 - 1:2000

### Immunogen Information

**Gene ID**

90957

**Swiss Prot**

Q6P158

**Immunogen**

Recombinant fusion protein containing a sequence corresponding to amino acids 700-900 of human DHX57 (NP\_945314.1).

**Synonyms**

DDX57; DHX57

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

**Source**

Rabbit

**Isotype**

IgG

**Purification**


Affinity purification

**Storage**

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of various lysates using DHX57 Rabbit pAb (A16589) at 1:1000 dilution.  
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.  
Lysates/proteins: 25μg per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.  
Detection: ECL Basic Kit (RM00020).  
Exposure time: 90s.