

A18666

Leader in Biomolecular Solutions for Life Science



YEATS2 Rabbit pAb

Catalog No.: A18666

1 Publications

Basic Information

Observed MW

151kDa

Calculated MW

151kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Mouse, Rat

Background

Summary: The protein encoded by this gene is a scaffolding subunit of the ATAC complex, which is a complex with acetyltransferase activity on histones H3 and H4. Alternative splicing results in multiple transcript variants encoding different isoforms.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

55689

Swiss Prot

Q9ULM3

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 130-331 of human YEATS2 (NP_060493.3).

Synonyms

FAME4; YEATS2

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.

Validation Data



Western blot analysis of various lysates using YEATS2 Rabbit pAb (A18666) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Enhanced Kit (RM00021).
Exposure time: 3min.