

A20228

Leader in Biomolecular Solutions for Life Science



CD31/PECAM1 Rabbit mAb

Catalog No.: A20228 **Recombinant** **1 Publications**

Basic Information

Observed MW

130-140kDa

Calculated MW

83KDa

Category

SMab Recombinant Monoclonal Antibody

Applications

WB,IHC-P,IP,ELISA

Cross-Reactivity

Human

CloneNo number

ARC5013-19

Background

The protein encoded by this gene is found on the surface of platelets, monocytes, neutrophils, and some types of T-cells, and makes up a large portion of endothelial cell intercellular junctions. The encoded protein is a member of the immunoglobulin superfamily and is likely involved in leukocyte migration, angiogenesis, and integrin activation.

Recommended Dilutions

WB 1:500 - 1:1000

IHC-P 1:50 - 1:200

IP 0.5µg-4µg antibody for
200µg-400µg extracts of
whole cells

Immunogen Information

Gene ID

5175

Swiss Prot

Q08481

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 630-738 of CD31/PECAM1 (NP_000433.4).

Synonyms

CD31; PECA1; GPIIA¹; PECAM-1; endoCAM; CD31/EndoCAM; CD31/PECAM1

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage


Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal,0.05% BSA,50% glycerol,pH7.3.


Validation Data




Western blot analysis of various lysates using CD31/PECAM1 Rabbit mAb (A20228) at 1:1000 dilution.
 Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
 Lysates/proteins: 25µg per lane.
 Blocking buffer: 3% nonfat dry milk in TBST.
 Detection: ECL Basic Kit (RM00020).
 Exposure time: 3s.



Immunohistochemistry analysis of CD31/PECAM1 in paraffin-embedded human kidney using CD31/PECAM1 Rabbit mAb (A20228) at dilution of 1:200 (40x lens). Perform high pressure antigen retrieval with 10 mM Tris/EDTA buffer pH 9.0 before commencing with IHC staining protocol.



Immunohistochemistry analysis of CD31/PECAM1 in paraffin-embedded human tonsil using CD31/PECAM1 Rabbit mAb (A20228) at dilution of 1:200 (40x lens). Perform high pressure antigen retrieval with 10 mM Tris/EDTA buffer pH 9.0 before commencing with IHC staining protocol.



Immunoprecipitation analysis of 300 µg extracts of Jurkat cells using 3 µg CD31/PECAM1 antibody (A20228). Western blot was performed from the immunoprecipitate using CD31/PECAM1 antibody (A20228) at a dilution of 1:1000.