

A2375

Leader in Biomolecular Solutions for Life Science



## Asymmetric DiMethyl-Histone H3-R26 Rabbit pAb

Catalog No.: A2375

1 Publications

### Basic Information

#### Observed MW

14kDa

#### Calculated MW

15kDa

#### Category

Polyclonal Antibody

#### Applications

WB,IF/ICC,ELISA

#### Cross-Reactivity

Human,Mouse,Rat,Other (Wide Range Predicted)

### Background

Histones are basic nuclear proteins that are responsible for the nucleosome structure of the chromosomal fiber in eukaryotes. Two molecules of each of the four core histones (H2A, H2B, H3, and H4) form an octamer, around which approximately 146 bp of DNA is wrapped in repeating units, called nucleosomes. The linker histone, H1, interacts with linker DNA between nucleosomes and functions in the compaction of chromatin into higher order structures. This gene is intronless and encodes a replication-dependent histone that is a member of the histone H3 family. Transcripts from this gene lack polyA tails but instead contain a palindromic termination element. This gene is found in the small histone gene cluster on chromosome 6p22-p21.3.

### Recommended Dilutions

|        |                |
|--------|----------------|
| WB     | 1:500 - 1:1000 |
| IF/ICC | 1:50 - 1:200   |

### Immunogen Information

**Gene ID**  
8290/8350

**Swiss Prot**  
Q16695/P68431

#### Immunogen

A synthetic asymmetric dimethylated peptide around R26 of human histone H3 (NP\_003520.1).

#### Synonyms

H3/j; H3C1; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FJ; H3C10; H3C11; HIST1H3J; Asymmetric DiMethyl-Histone H3-R26

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

**Source**  
Rabbit


**Isotype**  
IgG

**Purification**  
Affinity purification


#### Storage

Store at -20°C. Avoid freeze / thaw cycles.  
Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.

## Validation Data



Western blot analysis of various lysates using Asymmetric DiMethyl-Histone H3-R26 Rabbit pAb (A2375) at 1:1000 dilution.  
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.  
Lysates/proteins: 25µg per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.  
Detection: ECL Basic Kit (RM00020).  
Exposure time: 150s.



Immunofluorescence analysis of 293T cells using Asymmetric DiMethyl-Histone H3-R26 Rabbit pAb (A2375). Blue: DAPI for nuclear staining.