

Pan Acetyl-Lysine pAb

Catalog No	A2391	Category	Acetylated Antibodies
Applications	WB,IP	Observed MW	Refer to Figures
Cross-Reactivity	ALL	Calculated MW	

Immunogen Information

Immunogen	Acetylated BSA
Gene ID	
Swiss Prot	
Synonyms	


Product information

Source	Rabbit
Isotype	IgG
Purity	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB	1:500 - 1:1000
IP	1:20 - 1:50



Western blot analysis of extracts of HeLa cells, using Pan Acetyl-Lysine antibody (A2391).
 Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
 Lysates/proteins: 25ug per lane.
 Blocking buffer: 3% nonfat dry milk in TBST.