

A6442

Leader in Biomolecular Solutions for Life Science



UBE3C Rabbit pAb

Catalog No.: A6442

2 Publications

Basic Information

Observed MW

124kDa

Calculated MW

124kDa

Category

Mouse Monoclonal Antibody

Applications

WB,IF/ICC,ELISA

Cross-Reactivity

Human

Background

Enables ubiquitin protein ligase activity. Involved in protein polyubiquitination. Predicted to be located in nucleus. Predicted to be part of proteasome complex.

Recommended Dilutions

WB	1:500 - 1:2000
IF/ICC	1:10 - 1:100

Immunogen Information

Gene ID

9690

Swiss Prot

Q15386

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 970-1083 of human UBE3C (NP_055486.2).

Synonyms

RAUL; HECTH2; NDSMBA; NEDSMBA; UBE3C

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.