

A7530

Leader in Biomolecular Solutions for Life Science



DHRS3 Rabbit pAb

Catalog No.: A7530

1 Publications

Basic Information

Observed MW

33kDa

Calculated MW

34kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse, Rat

Background

Predicted to enable NAD-retinol dehydrogenase activity. Predicted to be involved in regulation of retinoic acid receptor signaling pathway and retinoid metabolic process. Predicted to act upstream of or within several processes, including animal organ morphogenesis; negative regulation of retinoic acid receptor signaling pathway; and regulation of ossification. Predicted to be located in endoplasmic reticulum membrane and photoreceptor outer segment membrane. Predicted to be active in lipid droplet.

Recommended Dilutions

WB 1:1000 - 1:2000

Immunogen Information

Gene ID

9249

Swiss Prot

O75911

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 30-150 of human DHRS3 (NP_004744.2).

Synonyms

SDR1; RDH17; Rsdr1; DD83.1; SDR16C1; retSDR1; DHRS3

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.

Validation Data



Western blot analysis of various lysates using DHRS3 Rabbit pAb (A7530) at 1:2000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 30s.