

AC028

Leader in Biomolecular Solutions for Life Science



HRP-conjugated β -Actin Rabbit mAb

Catalog No.: AC028

Recombinant

46 Publications

Basic Information

Observed MW

42kDa

Calculated MW

42kDa

Category

SMab Recombinant Monoclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse, Rat, Chicken, Zebrafish, Pig

CloneNo number

ARC5115-01

Conjugate

HRP

Background

This gene encodes one of six different actin proteins. Actins are highly conserved proteins that are involved in cell motility, structure, integrity, and intercellular signaling. The encoded protein is a major constituent of the contractile apparatus and one of the two nonmuscle cytoskeletal actins that are ubiquitously expressed. Mutations in this gene cause Baraitser-Winter syndrome 1, which is characterized by intellectual disability with a distinctive facial appearance in human patients. Numerous pseudogenes of this gene have been identified throughout the human genome.

Recommended Dilutions

WB 1:2000 - 1:50000

Immunogen Information

Gene ID

60

Swiss Prot

P60709

Immunogen

Recombinant protein of human HRP-conjugated β -Actin

Synonyms

BRWS1; PS1TP5BP1; HRP-conjugated β -Actin

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage


Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.05% proclin300, 0.05% BSA, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of various lysates using HRP-conjugated β-Actin Rabbit mAb (AC028) at 1:50000 dilution.
Lysates/proteins: 25μg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 3s.



Western blot analysis of lysates from HeLa, using HRP-conjugated β-Actin Rabbit mAb (AC028) at 1:1000-1:512000 dilution.
Lysates/proteins: 25μg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 1s.