

Rabbit anti Myc-Tag pAb

Catalog No	AE009	Category	Tag Antibodies
Applications	WB	Observed MW	Refer to Figures
Cross-Reactivity		Calculated MW	

Immunogen Information

Immunogen Recombinant protein of MYC-tag

Gene ID

Swiss Prot

Synonyms

Product information

Source Rabbit

Isotype IgG

Purity Affinity purification

Storage Store at -20°C. Avoid freeze / thaw cycles.
Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:1000 - 1:10000