

AP0778

Leader in Biomolecular Solutions for Life Science



## Phospho-SHIP1-Y1020 Rabbit pAb

Catalog No.: AP0778

### Basic Information

**Observed MW**

145kDa

**Calculated MW**

133kDa

**Category**

Polyclonal Antibody

**Applications**

WB, ELISA

**Cross-Reactivity**

Human

### Background

This gene is a member of the inositol polyphosphate-5-phosphatase (INPP5) family and encodes a protein with an N-terminal SH2 domain, an inositol phosphatase domain, and two C-terminal protein interaction domains. Expression of this protein is restricted to hematopoietic cells where its movement from the cytosol to the plasma membrane is mediated by tyrosine phosphorylation. At the plasma membrane, the protein hydrolyzes the 5' phosphate from phosphatidylinositol (3,4,5)-trisphosphate and inositol-1,3,4,5-tetrakisphosphate, thereby affecting multiple signaling pathways. The protein is also partly localized to the nucleus, where it may be involved in nuclear inositol phosphate signaling processes. Overall, the protein functions as a negative regulator of myeloid cell proliferation and survival. Mutations in this gene are associated with defects and cancers of the immune system. Deficiencies in the encoded protein, SHIP1, have been associated with Inflammatory Bowel Disease types such as Crohn's Disease and Ulcerative Colitis. Alternative splicing of this gene results in multiple transcript variants.

### Recommended Dilutions

WB 1:500 - 1:1000

### Immunogen Information

**Gene ID**

3635

**Swiss Prot**

Q92835

**Immunogen**

A synthetic phosphorylated peptide around Y1020 of human SHIP1 (NP\_001017915.1).

**Synonyms**

SHIP; SHIP1; SHIP-1; hp51CN; SIP-145; p150Ship; Phospho-SHIP1-Y1020

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

**Source**

Rabbit

**Isotype**

IgG

**Purification**


Affinity purification

**Storage**

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.

## Validation Data



Western blot analysis of lysates from Raji cells, using Phospho-SHIP1-Y1020 Rabbit pAb (AP0778) at 1:1000 dilution. Raji cells were treated by Pervanadate (1 mM) at 37°C for 30 minutes.  
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.  
Lysates/proteins: 25µg per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.  
Detection: ECL Basic Kit (RM00020).  
Exposure time: 1s.