

FITC Goat Anti-Mouse IgG (H+L)

| | | | |
|-------------------------|--------|----------------------|----------------------|
| Catalog No | AS001 | Category | Secondary Antibodies |
| Applications | ICC,FC | Observed MW | |
| Cross-Reactivity | | Calculated MW | |

Immunogen Information

| | |
|-------------------|-----------|
| Immunogen | Mouse IgG |
| Gene ID | |
| Swiss Prot | |
| Synonyms | |

Product information

| | |
|----------------|--|
| Source | Goat |
| Isotype | Fluorescein conjugated IgG |
| Purity | Affinity purification |
| Storage | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3. |

Background

Recommended Dilutions

| | |
|-----|--------------|
| ICC | 1:20 - 1:100 |
| FC | 1:20 - 1:100 |