

HRP Goat Anti-Mouse IgG (H+L)

| | | | |
|-------------------------|-----------|----------------------|----------------------|
| Catalog No | AS003 | Category | Secondary Antibodies |
| Applications | Elisa, WB | Observed MW | |
| Cross-Reactivity | | Calculated MW | |

Immunogen Information

Immunogen Mouse IgG

Gene ID

Swiss Prot

Synonyms

Product information

Source Goat

Isotype Horseradish peroxidase conjugated IgG

Purity Affinity purification

Storage Store at -20°C. Avoid freeze / thaw cycles.
Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

Elisa 1:5000 - 1:10000

WB 1:2000 - 1:10000