

AMCA Goat Anti-Rabbit IgG (H+L)

Catalog No	AS020	Category	Secondary Antibodies
Applications	ICC,FC	Observed MW	
Cross-Reactivity		Calculated MW	

Immunogen Information

Immunogen	Rabbit IgG
Gene ID	
Swiss Prot	
Synonyms	

Product information

Source	Goat
Isotype	AMCA conjugated IgG
Purity	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

ICC	1:50 - 1:200
FC	1:50 - 1:200