

HRP Goat Anti-Guinea Pig IgG (H+L)

Catalog No	AS025	Category	Secondary Antibodies
Applications	WB,IHC,ICC	Observed MW	
Cross-Reactivity		Calculated MW	

Immunogen Information

Immunogen	Guinea Pig IgG
Gene ID	
Swiss Prot	
Synonyms	

Product information

Source	Goat
Isotype	Horseradish peroxidase conjugated IgG
Purity	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB	1:2000 - 1:50000
IHC	1:200 - 1:2000
ICC	1:200 - 1:2000