

RP00056

Leader in Biomolecular Solutions for Life Science



Recombinant Human Sonic hedgehog protein N-product/SHH(C24IVI) Protein

Catalog No.: RP00056

Recombinant

Sequence Information

Species	Gene ID	Swiss Prot
Human	6469	Q15465

Tags

No tag

Synonyms

HHG1; HLP3; HPE3; MCOPCB5; SMMCI; TPT; TPTPS; SHH; HHG1; sonic hedgehog; HLP3; HPE3; MCOPCB5; SMMCI; TPT; TPTPS

Product Information

Source

<I>E. coli</I>

Purification

> 97% by SDS-PAGE.

Endotoxin

< 1.0 EU/μg of the protein by LAL method.

Formulation

Lyophilized from a 0.22 μm filtered solution of 20mM Tris, 300mM NaCl, pH 7.4. Contact us for customized product form or formulation.

Reconstitution

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Background

This protein is instrumental in patterning the early embryo. It has been implicated as the key inductive signal in patterning of the ventral neural tube, the anterior-posterior limb axis, and the ventral somites. Of three human proteins showing sequence and functional similarity to the sonic hedgehog protein of *Drosophila*, this protein is the most similar. The protein is made as a precursor that is autocatalytically cleaved; the N-terminal portion is soluble and contains the signalling activity while the C-terminal portion is involved in precursor processing. More importantly, the C-terminal product covalently attaches a cholesterol moiety to the N-terminal product, restricting the N-terminal product to the cell surface and preventing it from freely diffusing throughout the developing embryo. Defects in this protein or in its signalling pathway are a cause of holoprosencephaly (HPE), a disorder in which the developing forebrain fails to correctly separate into right and left hemispheres. HPE is manifested by facial deformities. It is also thought that mutations in this gene or in its signalling pathway may be responsible for VACTERL syndrome, which is characterized by vertebral defects, anal atresia, tracheoesophageal fistula with esophageal atresia, radial and renal dysplasia, cardiac anomalies, and limb abnormalities.

Basic Information

Description

Recombinant Human Sonic hedgehog protein N-product/SHH(C24IVI) Protein is produced by <I>E. coli</I> expression system. The target protein is expressed with sequence (Cys24-Gly197 (Cys24Ile-Val-Ile)) of human Sonic hedgehog N-product (Accession #NP_000184.1).

Bio-Activity

Measured by its ability to inhibit p53 expression in C3H10T1/2 mouse embryonic fibroblast cells. 1.25-2.5 μg/mL of Recombinant Human Sonic hedgehog can effectively decrease p53 expression.

Storage


Store the lyophilized protein at -20°C to -80 °C for long term.
After reconstitution, the protein solution is stable at -20 °C for 3 months, at 2-8 °C for up to 1 week. Avoid repeated freeze/thaw cycles.

Contact




www.abclonal.com


Validation Data



Recombinant Human Sonic hedgehog protein N-product/SHH(C24IVI) Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 22 kDa.



The purity of Human Sonic hedgehog N-product Protein (Cat.RP00056) was greater than 95% as determined by SEC-HPLC.



Recombinant Human Sonic hedgehog inhibits p53 expression in C3H10T1/2 mouse embryonic fibroblast cells. 1.25-2.5μg/mL of Recombinant Human Sonic hedgehog can effectively decrease p53 expression .