

RP00245

Leader in Biomolecular Solutions for Life Science



Recombinant Human Serum Amyloid P/SAP Protein

Catalog No.: RP00245

Recombinant

Sequence Information

Species	Gene ID	Swiss Prot
Human	325	P02743

Tags

C-His

Synonyms

APCS;HEL-S-92n;PTX2;SAP

Product Information

Source	Purification
HEK293 cells	> 97% by SDS-PAGE.

Endotoxin

< 0.01 EU/μg of the protein by LAL method.

Formulation

Lyophilized from a 0.22 μm filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation.

Reconstitution

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Contact



www.abclonal.com

Background

Serum amyloid P-component can interact with DNA and histones and may scavenge nuclear material released from damaged circulating cells. The identical serum form of amyloid P component (AP), a highly preserved plasma protein named for its ubiquitous presence in amyloid deposits. As a normal plasma protein first identified as the pentagonal constituent of in vivo pathological deposits called "amyloid". Serum amyloid P component represents another member of the pentraxin family, a highly conserved group of molecules that may play a role in innate immunity.

Basic Information

Description


Recombinant Human Serum Amyloid P/SAP Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (His 20 - Val 223) of human Serum Amyloid P (Accession #NP_001630.1) fused with a 6×His tag at the C-terminus.

Bio-Activity

Storage

Store the lyophilized protein at -20°C to -80 °C for long term. After reconstitution, the protein solution is stable at -20 °C for 3 months, at 2-8 °C for up to 1 week. Avoid repeated freeze/thaw cycles.

Validation Data



Recombinant Human Serum Amyloid P/SAP Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 26 kDa.