

RP00806

Leader in Biomolecular Solutions for Life Science



# Recombinant Mouse Fc gamma receptor IIB/CD32B Protein

Catalog No.: RP00806

Recombinant

## Sequence Information

Species	Gene ID	Swiss Prot
Mouse	14130	P08101

### Tags

C-6×His

### Synonyms

CD32;CD32B;FCG2;FCGR2;IGFR2;Fcgr2b

## Product Information

Source	Purification
HEK293 cells	> 95% by SDS-PAGE.

### Endotoxin

< 1 EU/μg of the protein by LAL method.

### Formulation

Lyophilized from a 0.2 μm filtered solution of PBS, pH7.4. Contact us for customized product form or formulation.

### Reconstitution

Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water.

## Background

Low affinity immunoglobulin gamma Fc region receptor II (CD32B) is a single-pass type I membrane protein and contains 2 Ig-like C2-type (immunoglobulin-like) domains. The inhibitory CD32B is expressed on B cells and myeloid dendritic cells. Ligation of CD32B on B cells downregulates antibody production and may, in some circumstances, promote apoptosis. Co-ligation of CD32B on dendritic cells inhibits maturation and blocks cell activation. CD32B may also be a target for monoclonal antibody therapy for malignancies.

## Basic Information

### Description

Recombinant Mouse Fc gamma receptor IIB/CD32B Protein is produced by Human cells expression system. The target protein is expressed with sequence (Thr30-Pro210) of mouse Fc gamma receptor IIB/CD32B (Accession #P08101) fused with a 6×His tag at the C-terminus.

### Bio-Activity

### Storage

Store the lyophilized protein at -20 °C to -80 °C for long term. After reconstitution, the protein solution is stable at -20 °C for 3 months, at 2-8 °C for up to 1 week. Avoid repeated freeze/thaw cycles.


## Contact



[www.abclonal.com](http://www.abclonal.com)

## Validation Data

---



Recombinant Mouse Fc gamma receptor IIB/CD32B Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.