

RP00981

Leader in Biomolecular Solutions for Life Science



Recombinant Human Semaphorin-4D/SEMA4D/CD100 Protein

Catalog No.: RP00981

Recombinant

Sequence Information

Species	Gene ID	Swiss Prot
Human	10507	Q92854

Tags

C-His

Synonyms

SEMA4D;C9orf164;CD100;M-sema-
G;SEMAJ;coll-4;Semaphorin 4D / CD100

Product Information

Source	Purification
HEK293 cells	> 97% by SDS- PAGE.

Endotoxin

< 0.1 EU/μg of the protein by LAL
method.

Formulation

Lyophilized from a 0.22 μm filtered
solution of PBS, pH 7.4. Contact us for
customized product form or formulation.

Reconstitution

Centrifuge the vial before opening.
Reconstitute to a concentration of
0.1-0.5 mg/mL in sterile distilled water.
Avoid vortex or vigorously pipetting the
protein. For long term storage, it is
recommended to add a carrier protein or
stabilizer (e.g. 0.1% BSA, 5% HSA, 10%
FBS or 5% Trehalose), and aliquot the
reconstituted protein solution to
minimize free-thaw cycles.

Background

This protein is a member of the semaphorin family of proteins and an important mediator of the movement and differentiation of multiple cell types, including those of the immune, vascular, and nervous systems. VEGF and SEMA4D had a positive correlation with the malignant degree of ovarian cancer, and SEMA4D can serve as an independent prognostic factor. SEMA4D was the first semaphorin described to have immune functions and serves important roles in T cell priming, antibody production, and cell-to-cell adhesion. Proteolytic cleavage of SEMA4D from the cell surface gives rise to a soluble fragment of SEMA4D (sSEMA4D). Similar to the transmembranal form, sSEMA4D is thought to have immunoregulatory properties.

Basic Information

Description

Recombinant Human Semaphorin-4D/SEMA4D/CD100 Protein is produced by HEK293 expression system. The target protein is expressed with sequence (Met22-Arg734) of human Semaphorin 4D (Accession #NP_006369.3) fused with a 6×His tag at the C-terminus.

Bio-Activity

Measured by its binding ability in a functional ELISA. Immobilized Human SEMA4D Protein at 1 μg/mL (100 μL/well) can bind CD100 Rabbit mAb with a linear range of 0.012-0.88 ng/mL.

Storage


Store the lyophilized protein at -20 °C to -80 °C for long term. After reconstitution, the protein solution is stable at -20 °C for 3 months, at 2-8 °C for up to 1 week. Avoid repeated freeze/thaw cycles.

Contact




www.abclonal.com

Validation Data



Recombinant Human Semaphorin-4D/SEMA4D/CD100 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 120 kDa.



Immobilized Human SEMA4D Protein at 1µg/mL (100 µL/well) can bind CD100 Rabbit mAb with a linear range of 0.012-0.88 ng/mL.