Leader in Biomolecular Solutions for Life Science

Recombinant SARS-CoV-2 Spike S1 Protein



Catalog No.: RP01262

Recombinant 7 Publications

Sequence Information

Gene ID Swiss Prot Species SARS-CoV-2 43740568

Tags

C-His

Synonyms

Envelope; SARS-CoV-2 Spike RBD (N501Y); Spike; Spike ECD; Spike RBD; Spike S1; Spike S2; Spike S2 ECD; S1-RBD protein; NCP-CoV RBD Protein; novel coronavirus RBD Protein;2019-nCoV RBD Protein; Sglycoprotein Subunit1 RBD Protein

Product Information

Source

Purification

HEK293 cells

>90% by SDS-PAGE;> 95% by

HPLC

Endotoxin

< 1.0 EU/µg of the protein by LAL method.

Formulation

Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4. or Supplied as a 0.22 µm filtered solution in PBS, pH 7.4. Contact us for customized product form or formulation.

Reconstitution

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid votex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stablizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Contact



www.abclonal.com

Background

The spike protein (S) of coronavirus (CoV) attaches the virus to its cellular receptor, angiotensin-converting enzyme 2 (ACE2). A defined receptor-binding domain (RBD) on S mediates thisinteraction. The S protein plays key parts in the induction of neutralizing-antibody and T-cellresponses, as well as protective immunity.

Basic Information

Description

Recombinant SARS-CoV-2(2019-nCoV) Spike S1 Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (Val16-Arg685) of SARS-COV-2(2019-nCoV) Spike S1 (Accession #YP 009724390.1) fused with a 6×His tag at the C-terminus.

Bio-Activity

1.Measured by its binding ability in a functional ELISA. Immobilized Recombinant SARS-CoV-2 Spike S1 at 2 µg/mL (100 µL/well) can bind Recombinant Human ACE2 with a linear range of 0.5-8.7 ng/mL.|2.Immobilized Human ACE2 on COOH Chip can bind SARS-COV-2 Spike S1 with an affinity constant of 11.4 nM as determined in a SPR assay (Nicoya OpenSPR).


Storage

Store the lyophilized protein at -20°C to -80 °C for long term.


-After reconstitution, the protein solution is stable at -20 °C for 3 months, at 2-8 °C for up to 1 week.

or This product is stable at \leq -70°C for up to 6 months from the date of receipt.
br/>For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze/thaw cycles.


Validation Data




Recombinant SARS-CoV-2 Spike S1 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 110-130 kDa.



than 95% as determined by SEC-HPLC.



Immobilized Recombinant SARS-COV-2 Spike S1 at 2 μ g/mL (100 μ L/well) can bind Recombinant Human ACE2 with a linear range of 0.5-8.7 ng/mL.



Immobilized Human ACE2 on COOH Chip, can bind SARS-COV-2 Spike S1 with an affinity constant of 11.4 nM as determined in a SPR assay (Nicoya OpenSPR).