

RP02410

Leader in Biomolecular Solutions for Life Science



# Biotinylated Recombinant Human IL-17A/CTLA-8 Protein

Catalog No.: RP02410

Recombinant

## Sequence Information

Species	Gene ID	Swiss Prot
Human	3605	Q16552-1

### Tags

C-His&Avi

### Synonyms

IL17A; CTLA-8; CTLA8; IL-17; IL-17A; IL17;  
interleukin-17A;CTLA-8;CTLA8;IL-17;IL-17 A;IL17

## Product Information

Source	Purification
HEK293 cells	> 95% by Tris-Bis PAGE

### Endotoxin

< 1 EU/μg of the protein by LAL method.

### Formulation

Lyophilized from a 0.22 μm filtered solution of PBS, pH 7.4.

### Reconstitution

Centrifuge the tube before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

## Background

## Basic Information

### Description

Biotinylated Recombinant Human IL17A Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Gly24-Ala155) of Human IL17A fused with a His and Avi tag at the C-terminal.

### Bio-Activity

### Storage

Store the lyophilized protein at -20°C to -80°C for long term. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week. Avoid repeated freeze/thaw cycles.


## Contact



[www.abclonal.com](http://www.abclonal.com)

## Validation Data

---



Biotinylated Human IL-17A Protein on Tris-Bis PAGE under reduced condition. The purity is greater than 95%.