

RP10009LQ

Leader in Biomolecular Solutions for Life Science



Recombinant Human linear Ub8 protein

Catalog No.: RP10009LQ

Recombinant

Sequence Information

Species	Gene ID	Swiss Prot
Human		

Tags
His Tag

Synonyms
linear Ub8

Product Information

Source	Purification
E. coli	~95% by SDS-PAGE

Endotoxin

Formulation
20 mM Tris, 150 mM NaCl, 2 mM β ME,
10% Glycerol

Reconstitution

Background

The linear ubiquitin chains are formed by conjugating the N-terminal methionine residue of one ubiquitin with the C-terminal glycine residue of another ubiquitin. This product is expressed in E. coli as a fusion protein.

Basic Information

Description

Bio-Activity

Storage


This product is stable at $\leq -70^{\circ}\text{C}$ for up to 6 months from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze/thaw cycles.

Contact



www.abclonal.com

Validation Data



Recombinant Human linear Ub8 protein was determined by SDS-PAGE with Coomassie Blue.