

RP10161LQ

Leader in Biomolecular Solutions for Life Science



# Recombinant Human Deubiquitinase OTUD6B protein

Catalog No.: RP10161LQ **Recombinant**

## Sequence Information

Species	Gene ID	Swiss Prot
Human	51633	Q8N6M0

### Tags

No Tag

### Synonyms

OTUD6B; CGI-77; DUBA-5; DUBA5; IDDFSD  
A

## Product Information

Source	Purification
<I>E. coli</I>	~ 90% by SDS-PAGE

### Endotoxin

### Formulation

20 mM Tris, 150 mM NaCl, 2 mM  $\beta$ ME,  
10% Glycerol

### Reconstitution

## Background

OTUD6B belongs to the OTU (ovarian tumour) family of deubiquitinating enzymes. Functions of this enzyme are not clear. In vitro, purified OTUD6B has not activity for cleavage of various di-Ub chains.

## Basic Information

### Description

### Bio-Activity

### Storage

This product is stable at  $\leq -70^{\circ}\text{C}$  for up to 6 months from the date of receipt. <br/>For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature.  
Avoid repeated freeze/thaw cycles.


## Contact



[www.abclonal.com](http://www.abclonal.com)

## Validation Data

---



Recombinant human Deubiquitinase  
OTUD6B protein was determined by SDS-  
PAGE with Coomassie Blue.