

RP10208K

Leader in Biomolecular Solutions for Life Science



Di-ubiquitin Explorer Kit

Catalog No.: RP10208K

Sequence Information

Species	Gene ID	Swiss Prot
Human		
proteins		

Tags

Synonyms

Product Information

Source	Purification
E. coli recombinant	≥95% by RP-HPLC and SDS-PAGE

Endotoxin

Formulation

20 mM Tris, 150 mM NaCl, 2 mM βME,
10% Glycerol, pH7.1 at 37 °C

Reconstitution

Background

This kit contains 5 µg each of linear, K6, K11, K27, K29, K33, K48, and K63-linked diubiquitin, a total of eight different diubiquitin forms. It can be used for determining linkage specificity of deubiquitinating enzymes and polyubiquitin-binding proteins

Basic Information

Description

Bio-Activity

Storage

This product is stable at ≤ -70°C for up to 6 months from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature.

Avoid repeated freeze/thaw cycles.

Contact



www.abclonal.com