

RP10230K

Leader in Biomolecular Solutions for Life Science



K63 Polyubiquitin Chain Capture Kit

Catalog No.: RP10230K

Sequence Information

| Species | Gene ID | Swiss Prot |
|----------|---------|------------|
| Human | | |
| proteins | | |

Tags

Synonyms

Product Information

| Source | Purification |
|---------------------|-------------------------|
| E. coli recombinant | Enzyme >95% by SDS-PAGE |

Endotoxin

Formulation

Enzyme stocked in 20 mM Tris, 150 mM NaCl, 2 mM β ME, 10% Glycerol

Reconstitution

Background

This kit is designed for enrichment of polyubiquitinated cellular proteins conjugated with lysine 63 using GST-TAB2 (NZF). Polyubiquitinated proteins bound on GST-TAB2 (NZF) can be precipitated using glutathione resin and then eluted by a buffer containing 10 mM glutathione. Polyubiquitinated proteins can be assayed by immunoblotting or mass spectrometry

Basic Information

Description

Bio-Activity

Storage

This product is stable at $\leq -70^{\circ}\text{C}$ for up to 6 months from the date of receipt.
For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature.

Avoid repeated freeze/thaw cycles.

Contact



www.abclonal.com