

# CYP1A2 Rabbit pAb

Catalog No.: A0062

## Basic Information

---

### Catalog No.

A0062

### Observed MW

54kDa

### Calculated MW

58kDa

### Category

Primary antibody

### Applications

WB

### Cross-Reactivity

Mouse, Rat

## Recommended Dilutions

---

WB 1:500 - 1:2000

## Background

---

This gene encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. The protein encoded by this gene localizes to the endoplasmic reticulum and its expression is induced by some polycyclic aromatic hydrocarbons (PAHs), some of which are found in cigarette smoke. The enzyme's endogenous substrate is unknown; however, it is able to metabolize some PAHs to carcinogenic intermediates. Other xenobiotic substrates for this enzyme include caffeine, aflatoxin B1, and acetaminophen. The transcript from this gene contains four Alu sequences flanked by direct repeats in the 3' untranslated region.

## Immunogen Information

---

### Gene ID

1544

### Swiss Prot

P05177

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 205-305 of human CYP1A2 (NP\_000752.2).

### Synonyms

CYP1A2;CP12;P3-450;P450(PA);CYP1A2

## Contact

---

 | [www.abclonal.com](http://www.abclonal.com)

## Product Information

---

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification


### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Validation Data

---



Western blot analysis of extracts of various cell lines, using CYP1A2 antibody (A0062) at 1:1000 dilution.

Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25ug per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 90s.