

DKK1 Rabbit pAb

Catalog No.: A0567 **1 Publications**

Basic Information

Catalog No.

A0567

Observed MW

36kDa

Calculated MW

28kDa

Category

Primary antibody

Applications

WB, IHC

Cross-Reactivity

Human, Mouse

Recommended Dilutions

WB	1:200 - 1:500
IHC	1:50 - 1:100

Background

This gene encodes a protein that is a member of the dickkopf family. It is a secreted protein with two cysteine rich regions and is involved in embryonic development through its inhibition of the WNT signaling pathway. Elevated levels of DKK1 in bone marrow plasma and peripheral blood is associated with the presence of osteolytic bone lesions in patients with multiple myeloma.

Immunogen Information

Gene ID	Swiss Prot
22943	O94907

Immunogen

Recombinant protein of human DKK1

Synonyms

DKK1;DKK-1;SK

Contact

 | www.abclonal.com

Product Information

Source	Isotype	Purification
Rabbit	IgG	Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.