

EIF3A Rabbit pAb

Catalog No.: A0573 **1 Publications**

Basic Information

Catalog No.

A0573

Observed MW

180kDa

Calculated MW

162kDa/166kDa

Category

Primary antibody

Applications

WB, IHC, IF

Cross-Reactivity

Human, Mouse

Recommended Dilutions

WB	1:500 - 1:1000
IHC	1:50 - 1:100
IF	1:50 - 1:100

Contact

 | www.abclonal.com

Background

Immunogen Information

Gene ID

8661

Swiss Prot

Q14152

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human EIF3A (NP_003741.1).

Synonyms

EIF3A; EIF3; EIF3S10; P167; TIF32; eIF3-p170; eIF3-theta; p180; p185

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage


Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Validation Data



Western blot analysis of extracts of various cell lines, using EIF3A antibody (A0573) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (RM00020). Exposure time: 30s.



Immunohistochemistry of paraffin-embedded human stomach using EIF3A antibody (A0573) at dilution of 1:100 (40x lens).



Immunofluorescence analysis of U2OS cells using EIF3A antibody (A0573) at dilution of 1:100. Blue: DAPI for nuclear staining.