

# PRKAR1A Rabbit pAb

Catalog No.: A0906

## Basic Information

### Catalog No.

A0906

### Observed MW

48kDa

### Calculated MW

38kDa/42kDa

### Category

Primary antibody

### Applications

WB, IHC, IF

### Cross-Reactivity

Human, Mouse

## Recommended Dilutions

|            |                |
|------------|----------------|
| <b>WB</b>  | 1:500 - 1:2000 |
| <b>IHC</b> | 1:50 - 1:200   |
| <b>IF</b>  | 1:50 - 1:200   |

## Contact



[www.abclonal.com](http://www.abclonal.com)

## Background

cAMP is a signaling molecule important for a variety of cellular functions. cAMP exerts its effects by activating the cAMP-dependent protein kinase, which transduces the signal through phosphorylation of different target proteins. The inactive kinase holoenzyme is a tetramer composed of two regulatory and two catalytic subunits. cAMP causes the dissociation of the inactive holoenzyme into a dimer of regulatory subunits bound to four cAMP and two free monomeric catalytic subunits. Four different regulatory subunits and three catalytic subunits have been identified in humans. This gene encodes one of the regulatory subunits. This protein was found to be a tissue-specific extinguisher that down-regulates the expression of seven liver genes in hepatoma x fibroblast hybrids. Mutations in this gene cause Carney complex (CNC). This gene can fuse to the RET protooncogene by gene rearrangement and form the thyroid tumor-specific chimeric oncogene known as PTC2. A nonconventional nuclear localization sequence (NLS) has been found for this protein which suggests a role in DNA replication via the protein serving as a nuclear transport protein for the second subunit of the Replication Factor C (RFC40). Several alternatively spliced transcript variants encoding two different isoforms have been observed.

## Immunogen Information

### Gene ID

5573

### Swiss Prot

P10644

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human PRKAR1A (NP\_002725.1).

### Synonyms

PRKAR1A;ACRDYS1;ADOHR;CAR;CNC;CNC1;PKR1;PPNAD1;PRKAR1;TSE1

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification


Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Validation Data




Western blot analysis of extracts of various cell lines, using PRKAR1A antibody (A0906) at 1:1000 dilution.

Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25ug per lane.

Blocking buffer: 3% nonfat dry milk in TBST.



Western blot analysis of extracts from normal (control) and PRKAR1A knockout (KO) 293T cells, using PRKAR1A antibody (A0906) at 1:1000 dilution.


Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25ug per lane.


Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Enhanced Kit (RM00021).

Exposure time: 180s.



Immunohistochemistry of paraffin-embedded human oophoroma using PRKAR1A antibody (A0906) at dilution of 1:100 (40x lens).



Immunofluorescence analysis of NIH/3T3 cells using PRKAR1A antibody (A0906) at dilution of 1:100. Blue: DAPI for nuclear staining.