

# GNG7 Rabbit pAb

Catalog No.: A10009 **2 Publications**

## Basic Information

---

**Catalog No.**

A10009

**Observed MW**

12kDa

**Calculated MW**

7kDa

**Category**

Primary antibody

**Applications**

WB

**Cross-Reactivity**

Mouse, Rat

## Background

---

## Recommended Dilutions

---

**WB** 1:500 - 1:1000

## Immunogen Information

---

**Gene ID**

2788

**Swiss Prot**

O60262

**Immunogen**


Recombinant fusion protein containing a sequence corresponding to amino acids 1-68 of human GNG7 (NP\_443079.1).

**Synonyms**

GNG7

## Contact

---

 | [www.abclonal.com](http://www.abclonal.com)

## Product Information

---

**Source**

Rabbit

**Isotype**

IgG

**Purification**

Affinity purification


**Storage**

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Validation Data

---



Western blot analysis of extracts of various cell lines, using GNG7 antibody (A10009) at 1:1000 dilution.

Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25ug per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Enhanced Kit (RM00021).

Exposure time: 90s.