

# RTN4RL1 Rabbit pAb

Catalog No.: A10021

## Basic Information

---

### Catalog No.

A10021

### Observed MW

49kDa

### Calculated MW

49kDa

### Category

Primary antibody

### Applications

WB

### Cross-Reactivity

Mouse, Rat

## Recommended Dilutions

---

**WB** 1:1000 - 1:4000

## Background

---

## Immunogen Information

---

### Gene ID

146760

### Swiss Prot

Q86UN2

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 220-420 of human RTN4RL1 (NP\_848663.1).

### Synonyms

RTN4RL1;NGRH2;NgR3

## Contact

---

 | [www.abclonal.com](http://www.abclonal.com)

## Product Information

---

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification


### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

## Validation Data

---



Western blot analysis of extracts of various cell lines, using RTN4RL1 antibody (A10021) at 1:1000 dilution.

Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25ug per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 20s.