

LY9 Rabbit pAb

Catalog No.: A10032

Basic Information

Catalog No.

A10032

Observed MW

90kDa

Calculated MW

57kDa/70kDa/71kDa/72kDa

Category

Primary antibody

Applications

WB

Cross-Reactivity

Human, Mouse

Recommended Dilutions

WB 1:500 - 1:1000

Background

Immunogen Information

Gene ID

4063

Swiss Prot

Q9HBG7


Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 160-360 of human LY9 (NP_001248386.1).

Synonyms

LY9;CD229;SLAMF3;hly9;mLY9

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

Validation Data



Western blot analysis of extracts of various cell lines, using LY9 antibody (A10032) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit (RM00021). Exposure time: 90s.