

IRF9 Rabbit pAb

Catalog No.: A1004 **1 Publications**

Basic Information

Catalog No.

A1004

Observed MW

43kDa

Calculated MW

43kDa

Category

Primary antibody

Applications

WB

Cross-Reactivity

Human

Background

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

10379

Swiss Prot

Q00978

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 113-255 of human IRF9 (NP_006075.3).

Synonyms

IRF9;IRF-9;ISGF3;ISGF3G;p48

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

Validation Data



Western blot analysis of extracts of various cell lines, using IRF9 antibody (A1004) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25ug per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 60s.