

RPL8 Rabbit pAb

Catalog No.: A10042 **1 Publications**

Basic Information

Catalog No.

A10042

Observed MW

28kDa

Calculated MW

28kDa

Category

Primary antibody

Applications

WB, IHC, IF

Cross-Reactivity

Human, Mouse, Rat

Recommended Dilutions

WB	1:1000 - 1:4000
IHC	1:50 - 1:200
IF	1:100 - 1:200

Contact

 | www.abclonal.com

Background

Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 60S subunit. The protein belongs to the L2P family of ribosomal proteins. It is located in the cytoplasm. In rat, the protein associates with the 5.8S rRNA, very likely participates in the binding of aminoacyl-tRNA, and is a constituent of the elongation factor 2-binding site at the ribosomal subunit interface. Alternatively spliced transcript variants encoding the same protein exist. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome.

Immunogen Information

Gene ID	Swiss Prot
6132	P62917

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-257 of human RPL8 (NP_000964.1).

Synonyms

RPL8;L8

Product Information


Source	Isotype	Purification
Rabbit	IgG	Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of extracts of various cell lines, using RPL8 antibody (A10042) at 1:1000 dilution.


Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25ug per lane.


Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).


Exposure time: 5s.




Immunohistochemistry of paraffin-embedded rat brain using RPL8 Rabbit pAb (A10042) at dilution of 1:50 (40x lens).




Immunohistochemistry of paraffin-embedded mouse brain using RPL8 Rabbit pAb (A10042) at dilution of 1:50 (40x lens).



Immunohistochemistry of paraffin-embedded mouse heart using RPL8 Rabbit pAb (A10042) at dilution of 1:50 (40x lens).



Immunofluorescence analysis of HeLa cells using RPL8 antibody (A10042) at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.



Confocal immunofluorescence analysis of HeLa cells using RPL8 antibody (A10042) at dilution of 1:200. Blue: DAPI for nuclear staining.