

CEACAM4 Rabbit pAb

Catalog No.: A10055

Basic Information

Catalog No.

A10055

Observed MW

24kDa

Calculated MW

25kDa

Category

Primary antibody

Applications

WB

Cross-Reactivity

Human

Recommended Dilutions

WB 1:500 - 1:1000

Background

Immunogen Information

Gene ID

1089

Swiss Prot

O75871

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 36-155 of human CEACAM4 (NP_001808.2).

Synonyms

CEACAM4;CGM7;CGM7_HUMAN;NCA

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

Validation Data



Western blot analysis of extracts of various cell lines, using CEACAM4 antibody (A10055) at 1:1000 dilution.

Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25ug per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Enhanced Kit (RM00021).

Exposure time: 20s.