

FCGR2C Rabbit pAb

Catalog No.: A10062

Basic Information

Catalog No.

A10062

Observed MW

45kDa

Calculated MW

25kDa/29kDa/30kDa/35kDa

Category

Primary antibody

Applications

WB

Cross-Reactivity

Human

Recommended Dilutions

WB 1:1000 - 1:4000

Background

Immunogen Information

Gene ID

9103

Swiss Prot

P31995

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 43-223 of human FCGR2C (NP_963857.3).

Synonyms

FCGR2C;CD32;CD32C;CDW32;FCG2;FCRIIC;IGFR2

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

Validation Data



Western blot analysis of extracts of various cell lines, using FCGR2C antibody (A10062) at 1:1000 dilution.

Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25ug per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 10s.