

LY75 Rabbit pAb

Catalog No.: A10065

Basic Information

Catalog No.

A10065

Observed MW

205kDa

Calculated MW

19kDa/26kDa/198kDa/209kDa/215kDa

Category

Primary antibody

Applications

WB

Cross-Reactivity

Mouse

Recommended Dilutions

WB 1:500 - 1:1000

Background

Immunogen Information

Gene ID

4065

Swiss Prot

O60449

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 440-680 of human LY75 (NP_002340.2).

Synonyms

LY75;CD205;CLEC13B;DEC-205;GP200-MR6;LY-75

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

Validation Data



Western blot analysis of extracts of mouse lung, using LY75 antibody (A10065) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25ug per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 90s.