

INSRR Rabbit pAb

Catalog No.: A10083

Basic Information

Catalog No.

A10083

Observed MW

105kDa

Calculated MW

143kDa

Category

Primary antibody

Applications

WB

Cross-Reactivity

Mouse, Rat

Recommended Dilutions

WB 1:500 - 1:1000

Background

Immunogen Information

Gene ID

3645

Swiss Prot

P14616

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 747-921 of human INSRR (NP_055030.1).

Synonyms

INSRR;IRR

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of extracts of various cell lines, using INSRR antibody (A10083) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit (RM00021). Exposure time: 90s.