#### A10098

Leader in Biomolecular Solutions for Life Science

# **COX IV Mouse mAb**

Catalog No.: A10098 6 Publications



### Background

Cytochrome c oxidase (COX) is the terminal enzyme of the mitochondrial respiratory chain. It is a multi-subunit enzyme complex that couples the transfer of electrons from cytochrome c to molecular oxygen and contributes to a proton electrochemical gradient across the inner mitochondrial membrane. The complex consists of 13 mitochondrial- and nuclear-encoded subunits. The mitochondrially-encoded subunits perform the electron transfer and proton pumping activities. The functions of the nuclear-encoded subunits are unknown but they may play a role in the regulation and assembly of the complex. This gene encodes the nuclear-encoded subunit IV isoform 1 of the human mitochondrial respiratory chain enzyme. It is located at the 3' of the NOC4 (neighbor of COX4) gene in a head-to-head orientation, and shares a promoter with it. Pseudogenes related to this gene are located on chromosomes 13 and 14. Alternative splicing results in multiple transcript variants encoding different isoforms.

Swiss Prot P13073

# 17kDa

**Basic Information** 

Calculated MW 20kDa

**Observed MW** 

Category SMab Recombinant Monoclonal Antibody

Applications WB,IHC-P,IF/ICC,IP,FC,ELISA

Cross-Reactivity Human,Mouse,Rat,Hamster,Monkey,Goat

## **Recommended Dilutions**

	mmun	ogen	Infor	mation
-				

WB	1:1000 - 1:4000
ІНС-Р	1:50 - 1:200
IF/ICC	1:50 - 1:200
IP	1:20 - 1:50
FC	1:20 - 1:50

1327		
Immunogen		

A synthetic peptide of human COX IV

Synonyms

COX4; COXIV; COX4-1; COXIV-1; MC4DN16; COX IV-1; COX IV

#### Contact

### **Product Information**

	www.abclonal.com	Source
•		Mouse


**Isotype** IgG Purification Affinity purification

Storage

Gene ID

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

## Validation Data



Western blot analysis of various lysates using COX IV Mouse mAb (A10098) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Mouse IgG (H+L) (AS003) at 1:10000 dilution. Lysates/proteins: 25µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (RM00020). Exposure time: 30s.